



STAINLESS STEEL SANITARY HOUSING SINGLE CARTRIDGE HOUSING

FEATURES:

- Material: 316L Stainless Steel
- Operating Pressure: 150 PSI
- Operating Temperature: 80°C (175°F) Max
- Finish: Electropolished, 25 RA External Finish, 20 RA Internal Finish
- Inlet/Outlet Type: Tri-clover
- Cartridge Connector:
 - 222 Style
 - 226 Style
- Model: NF Single Cartridge Stainless Steel Sanitary Housing



BUILDING A PART NUMBER

Model	# of Cartridges		Length	Cartridge Connectors
NF	1	R	30	
	1 = Single Round		10 = 10" 20 = 20" 30 = 30"	222 Style = 222 226 Style = BLANK

*Custom materials or sizes may be available upon request